



PACIFIC SHOWER BATHING TROLLEY



SKU: 3100002

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Unique solution for safe patient showering, bathing and transfer. The versatile Pacific Shower-Bathing Trolley provides a safe environment for showering, bathing and transfer.

The soft, padded top is ideal for patients with delicate skin or disabilities.

The electric high-low and unique electric tub-side actuation provides an ergonomic working height for caregivers and ensures safer patient transfers.

The Pacific Shower-Bathing Trolley has been carefully designed to be intuitive and easy to use meaning minimal user training is required.

The one piece seamless design also makes it easy to clean.

FEATURES

- Fully electric controls with battery backup
- Tub sides lower to become flush with patient platform
- High performance LINAK and TENTE components
- 200kg Safe Working Load
- 125mm corrosion resistant castors
- Safe 24V DC system with remote battery charger
- Soft, padded polyurethane tub
- Low battery warning light on handset
- Waterproof pillow
- Push handles at both ends
- Siderail (optional)



SPECIFICATIONS

Maximum SWL: 200kg Trolley Mass: 100kg

Overall Dimensions: 2165 x 775mm

Tub Platform Dimension: 1825 x 480mm Wheel Base: 1165mm, wheel track; 690mm Castor Diameter: 125mm single wheel Maximum Tub Depth (drain end): 210mm Maximum Tub Depth (head end): 180mm

Power Source: Rechargeable battery, 24V DC, 2.9Ah

Battery Charger Input: 110 or 240V AC, 50-60Hz, 650mA (max)

Battery Charge Time: 6 hours

Number of Full Load Cycles: Approx. 30

Time To Raise or Lower Patient Platform: 30 seconds

Time to raise or lower Tub Sides: 22 seconds

CLEANING

<u>View the approved cleaning products for mattress covers, upholstery tops, and shower trolley tub</u> fabrics here.

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 319903 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au



NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au