

SPIRIT MEDICAL PENLIGHT BI-FUNCTIONAL

\$62.00



SKU: SPT000



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

This model of diagnostic penlight torch offers a sturdy lightweight design, as well as a two way on / off switch for convenient use. The "TipTop" push button switch allows the unit to be activated in short bursts, while the pocket clip can lock the unit in the 'on' position. This allows both convenience of use as well as increased battery life.

FEATURES

- Two x AAA batteries included
- LED light
- Lightweight
- Fitted with LED light bulb for enhanced light quality this is an outstanding medical penlight.
- Suitable for optical use