



## **PERFORMANCE HEALTH KING SHOWER CHAIR WITH CUTOUT, BACK & ARMRESTS**



**SKU:** KING-BSC1BA-01KN-CBR4-CSC4

## **PRODUCT DESCRIPTION**

### **Purpose and benefits**

The King shower chair with cutout provides comfort, safety, and practical support during showering. It is lightweight and easy to adjust in height, allowing a suitable fit for different users. Non slip rubber feet keep the chair steady on wet or slippery surfaces.

The chair is made from durable plastic with sturdy aluminium legs for long term use. The seat features an anti slip surface and a secure locking system. A wide seat with smooth edges adds comfort during longer sitting periods.

The U shaped cutout can be placed at the front or back to support personal care needs. Backrest and armrests provide additional stability and comfort. The chair arrives flat packed for convenient shipping and storage, assembles quickly without tools, and all components are dishwasher safe for hygienic cleaning.

### **Features**

- Lightweight design with height adjustability
- Non slip rubber feet for secure positioning
- Durable plastic seat with sturdy aluminium legs
- Wide seat with smooth, rounded edges
- Moveable U shaped cutout for flexible positioning
- Backrest and armrests for added support
- Tool free assembly with dishwasher safe parts
- Shower chair and cushion sold separately

## Specifications

Max User Weight	150 kg
Seat Width	46 cm
Seat Depth	32.5 cm
Seat Height	39 to 54 cm
Width Between Armrests	45.5 cm
Armrest Height	23 cm
Backrest Height	36 cm
Backrest Length	39 cm
Total Width	53.5 cm
Total Depth	42 cm
Certifications	Manufactured with biocompatible materials according to EN ISO 10993 5 2009 and EN ISO 10993 10 2013

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)