



PERFORMANCE HEALTH SAVANAH SLATTED SHOWER BOARD



SKU: PAT-AA1091B-PARENT

PRODUCT DESCRIPTION

Purpose & Benefits

The **Savanah Slatted Shower Board** is designed to provide a stable seating surface for users who shower over the bath. Featuring two additional slats compared to the Savanah Slatted Bath Board, it offers increased comfort, support, and a larger surface area for greater confidence during transfers. The board is suitable for wider baths and features adjustable fittings to accommodate a range of bath widths. Fully sealed slats and stainless steel fittings provide corrosion resistance, making it suitable for regular use in wet environments. It is also compatible with the Savanah Slatted Bath Board Handle for additional support.

Features

- Designed for showering over the bath
- Two additional slats provide increased comfort and support
- Large seating surface improves confidence during transfers
- Suitable for wider baths
- Corrosion-resistant with fully sealed slats and stainless steel fittings
- Compatible with the Savanah Slatted Bath Board Handle

- Maximum user weight of 190 kg

Specifications

Size	711 mm
Board Dimensions	71.1 × 23.2 cm
Adjustable Bath Width	53.3 cm to 63.5 cm
Size	762 mm
Board Dimensions	76.2 × 23.2 cm
Adjustable Bath Width	58.4 cm to 68.6 cm
Maximum User Weight	190 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au

