


## PERMOBIL ROHO BARIATRIC MATTRESS SYSTEM WITH COVER



**SKU:** R-MATT-BARI42SYSC

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Permobil ROHO Bariatric Mattress System** with cover is a non-powered, air-filled pressure care solution engineered to deliver exceptional pressure redistribution and skin protection for bariatric users. Utilising ROHO's DRY FLOATATION® technology, it conforms precisely to each individual's body shape, reducing the risk of pressure injuries while promoting wound healing and comfort. Its modular and adaptable design makes it suitable for both clinical and home care environments where bariatric pressure management is a priority.

### Features

- Advanced DRY FLOATATION® technology with interconnected air cells for superior pressure redistribution
- Non-powered, air-filled design eliminates noise, heat, and electrical dependence
- Dynamic air cells adjust individually to conform to the user's body contours and movement
- Modular configuration allows width expansion to accommodate bariatric bed frames
- Approximately 8-9 cm (3.5") high air cells for optimal immersion and stability
- Durable neoprene rubber construction – flame resistant, easy to clean, and long lasting

- Drainage holes assist with moisture control for a hygienic and comfortable surface
- Includes washable, protective cover to maintain hygiene and prolong product life

## Specifications

Specification	Details
Model	Permobil ROHO Bariatric Mattress System
Technology	ROHO DRY FLOATATION® air cell pressure redistribution
Power Source	Non-powered (silent operation)
Cell Height	Approximately 8–9 cm (3.5")
Material	Flame-resistant neoprene rubber
Cover	Washable, removable, and breathable
Drainage	Integrated drainage holes to remove moisture
Customisation	Modular width expansion compatible with standard ROHO overlays
User Type	Designed for bariatric patients requiring pressure injury prevention or wound care

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)