



PERMOBIL ROHO BARIATRIC MATTRESS SYSTEM WITH COVER



SKU: R-MATT-BARI42SYSC

PRODUCT DESCRIPTION

Purpose & Benefits

The **Permobil ROHO Bariatric Mattress System** with cover is a non-powered, air-filled pressure care solution engineered to deliver exceptional pressure redistribution and skin protection for bariatric users. Utilising ROHO's DRY FLOATATION® technology, it conforms precisely to each individual's body shape, reducing the risk of pressure injuries while promoting wound healing and comfort. Its modular and adaptable design makes it suitable for both clinical and home care environments where bariatric pressure management is a priority.

Features

- Advanced DRY FLOATATION® technology with interconnected air cells for superior pressure redistribution
- Non-powered, air-filled design eliminates noise, heat, and electrical dependence
- Dynamic air cells adjust individually to conform to the user's body contours and movement
- Modular configuration allows width expansion to accommodate bariatric bed frames
- Approximately 8-9 cm (3.5") high air cells for optimal immersion and stability
- Durable neoprene rubber construction flame resistant, easy to clean, and long lasting



- Drainage holes assist with moisture control for a hygienic and comfortable surface
- Includes washable, protective cover to maintain hygiene and prolong product life

Specifications

Specification Details

Model Permobil ROHO Bariatric Mattress System

Technology ROHO DRY FLOATATION® air cell pressure redistribution

Power Source Non-powered (silent operation) Approximately 8-9 cm (3.5") Cell Height

Material Flame-resistant neoprene rubber

Cover Washable, removable, and breathable

Integrated drainage holes to remove moisture Drainage

Customisation Modular width expansion compatible with standard ROHO overlays

Designed for bariatric patients requiring pressure injury prevention or wound User Type

care