

## **PETA MULTI-PURPOSE FOOD PREPARATION BOARD**

**\$199.00**



**SKU: PKT-FPB**



## **PRODUCT DESCRIPTION**

### **Purpose & Benefits of Peta Multi-Purpose Food Preparation Board**

The **Peta Multi-Purpose Food Preparation Board** is a versatile kitchen aid designed to assist individuals with limited hand function, weakness, or the use of only one hand. Ideal for people living with arthritis, stroke recovery, or other conditions affecting grip and dexterity, this board makes food preparation safer, easier, and more independent.

Equipped with a range of features including stainless steel spikes, a clamp, and non-slip suction feet, the board holds items securely in place while cutting, peeling, grating, or spreading. Its durable and easy-to-clean design supports a variety of kitchen tasks, reducing the need for assistance during meal prep.

### **Key Features:**

- Built-in clamp holds bowls, jars, or food items in place
- Stainless steel spikes keep fruit, vegetables, or bread stable
- Non-slip suction feet for secure positioning on bench tops
- Ideal for one-handed use or users with limited hand strength
- Durable, easy-clean surface for everyday kitchen use

- Supports safe and independent food preparation
- Suitable for home kitchens, rehab centres, and aged care settings

**Specifications:**

- Length: 44 cm
- Width: 30 cm
- Height: 5.5 cm