




PLEXA HOMECARE BED



SKU: PLEXA-S-00000-NAT

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Plexa Series Bed is designed to provide safe and comfortable sleeping arrangements for individuals with complex needs. This bed includes full-length rails to provide security, featuring a locking system, ensuring the rails stay in place while the user is sleeping. An additional feature is the height adjustment making transfers in and out of bed with ease, designed with a durable frame and available in various sizes to accommodate the user's needs.

FEATURES

- APO Grey Powdercoated Frame
- Auto-regression Back Rest
- Trendelenburg
- 75mm Castors
- Vertical Lift
- Extremely Low Minimum Position
- Full-Length Timber Drop Down Rails
- Chair Position at the push of a button

SPECIFICATIONS

Height Range Max: 725mm – Min: 170mm

Weight Capacity 250kg

Angles

Back – 62 degrees

Knee – 35 degrees

Forward Tilt – 15 degrees

Reverse Tilt – 15 degrees

Single

Overall Dimensions 2150 x 1040mm

Sleeping Surface 2000 x 920mm

King Single

Overall Dimensions 2150 x 1210mm

Sleeping Surface 2000 x 1070mm

Double

Overall Dimensions 21650 x 1510mm

Sleeping Surface 2000 x 1370mm

BROCHURE

[Plexa Series Bed Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au