



PREMIUMLIFT DELUXE FLOOR BED



SKU: PLDFBK

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Deluxe Floorline Bed is the ultimate non-clinical fully functional floorline bed. It offers a unique Hidden Castors System that will revolutionise the home care market and truly enhances patient's recovery whilst promoting emotional comfort.

The advanced frame system can be lowered to just 18cm from the floor level and can elevate to an incredible height of 73cm and meet OH&S height standard requirements. Deluxe Floor Bed will promote the carer's wellbeing, prevent injury and offer a memorable patient experience and life changing care.

Uses & Benefits

Fall Prevention After Stroke Care

Alzheimer Care Injury Recovery

Carers Support Disability Care

FEATURES

- Min Height 18cm
- Hidden Castors System
- Non-Clinical Design
- Remote Control
- Five-function electric bed

- Zero Gravity & Seating Functions
- TGA Approved and ISO Certified
- 10 year limited warranty
- SWL 200kg

OPTIONS & ACCESSORIES

- Head & Foot Boards
- 20cm Length Extension
- Under Bed Light
- Foldable Side Rails
- Grab Rail
- Remote Control Pocket Holder
- Fall Mat
- Mattress
- External Battery

SPECIFICATIONS

King Single Size Mattress Platform	1070mm x 2050mm
Height Range	180mm – 730mm
Backrest Adjustament	0 – 70 Degrees
Legrest Adjustment	0 – 30 Degrees
Trendelenburg	20 Degrees
Safe Working Load	200 kg
Weight / Mass	80 kg
Locking Castors	50mm Twin wheel

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au