



PREMIUMLIFT DELUXE FLOOR BED



SKU: PLDFBK

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Deluxe Floorline Bed is the ultimate non-clinical fully functional floorline bed. It offers a unique Hidden Castors System that will revolutionise the home care market and truly enhances patient's recovery whilst promoting emotional comfort.

The advanced frame system can be lowered to just 18cm from the floor level and can elevate to an incredible height of 73cm and meet OH&S height standard requirements. Deluxe Floor Bed will promote the carer's wellbeing, prevent injury and offer a memorable patient experience and life changing care.

Uses & Benefits

Fall Prevention After Stroke Care

Alzheimer Care Injury Recovery

Carers Support Disability Care

FEATURES

- Min Height 18cm
- Hidden Castors System
- Non-Clinical Design
- Remote Control
- Five-function electric bed



- Zero Gravity & Seating Functions
- TGA Approved and ISO Certified
- 10 year limited warranty
- SWL 200kg

OPTIONS & ACCESSORIES

- Head & Foot Boards
- 20cm Length Extension
- Under Bed Light
- Foldable Side Rails
- Grab Rail
- Remote Control Pocket Holder
- Fall Mat
- Mattress
- External Battery

SPECIFICATIONS

King Single Size Mattress Platform

Height Range

Backrest Adjustament

Legrest Adjustment

Trendelenburg

Safe Working Load

Weight / Mass

Locking Castors

1070mm x 2050mm

180mm - 730mm

0 - 70 Degrees

0 - 30 Degrees

20 Degrees

200 kg

80 kg

50mm Twin wheel