


## PREMIUMLIFT ULTRA LOW BED



**SKU:** PLULBK

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

PremiumLift Ultra Low Bed is a unique system, which allows the platform to be lowered to as little as 6.7cm from the floor ensuring complete user and patient safety. PremiumLift Ultra Low Bed is the ultimate in floor level bed, designed in Germany to the higher quality specification and equipped with a state of the art Anti-Entrapment System.

### Uses & Benefits

Fall Prevention    After Stroke Care

Alzheimer Care    Injury Recovery

Carers Support    Disability Care

### FEATURES

- Min Height 6.7cm
- Anti-Entrapment System
- Five-function electric bed
- TGA Approved and ISO Certified
- 10 year limited warranty
- 4 section bed for easy assembly and dismantle
- 9V Battery Backup
- Beech Head & Foot Boards
- SWL 185kg

## OPTIONS AND ACCESSORIES

- King Single Width Extension Kit (1060mm)
- Foldable Side Rails
- Grab Rail
- Self Help Pole
- Fall Mat
- Mattress

**King Single Extension Kit** - Our unique King Single Width Extension Kit will allow the user to transform the Single size frame into a King Single size without having to buy a new bed. The extension kit is easily installed and removed.

## SPECIFICATIONS

<https://youtu.be/h4Z1bLWN548>