



## PRIDE LITESTREAM FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR



**SKU:** PWS654185

## PRODUCT DESCRIPTION

### Purpose & Benefits of Pride Litestream Folding Carbon Fibre Electric Wheelchair

The Litestream Folding Carbon Fibre Electric Wheelchair by Pride is designed for optimal mobility and convenience. This lightweight, versatile wheelchair is perfect for users seeking a travel-friendly mobility solution without compromising on performance. Its carbon fibre construction ensures durability and ease of handling, making it ideal for both indoor and outdoor use.

### Features

- **Carbon Fibre Frame:** Utilises high-strength yet lightweight carbon fibre materials to enhance manoeuvrability and transportability.
- **Compact Folding Design:** Easily folds to a compact size for convenient storage and transportation, suitable for car travel and air travel.
- **Advanced Drive System:** Features a powerful drive system capable of reaching speeds up to 6 km/h with a reliable range of 15 km on a single charge.
- **User Comfort:** Equipped with a padded seat and backrest, providing comfortable support for extended use.
- **Intuitive Controls:** Includes a user-friendly interface with a digital LCD that shows vital information such as speed, battery life, and travel distance.
- **Safety Features:** Comes with integrated safety elements, including an emergency stop function, anti-tip wheels, and reflective panels for increased visibility.

### Specifications

<b>Product Weight</b>	16 kg (ex. batteries)
<b>Safe Working Load</b>	136 kg
<b>Total Length</b>	930 mm
<b>Total Width</b>	577 mm
<b>Maximum Speed</b>	6 km/h
<b>Range</b>	Up to 15 km
<b>Seat width</b>	380 mm
<b>Batteries</b>	1 × 24V 12AH Lithium Battery (288wh)
<b>Turning Radius</b>	1200 mm

### Warranty:

- Frame & Electronic Components: 2 Years
- Batteries: 1 Year

The Litestream Folding Carbon Fibre Electric Wheelchair combines advanced technology with practical design elements to offer enhanced mobility and independence. Its robust construction and foldable nature make it a top choice for individuals looking for a portable and high-performance mobility solution.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)