

PRIDE PATHRIDER 140 XL SCOOTER



SKU: PR140XL-4 RED-H



PRODUCT DESCRIPTION

Purpose & Benefits of Pride Pathrider 140 XL Scooter

The **Pride Pathrider 140 XL Scooter** is a heavy-duty, outdoor-ready mobility scooter built for high performance, comfort, and reliability. Designed for users who require long-distance travel and rugged outdoor capability, the Pathrider 140 XL features a powerful motor, large pneumatic tyres, and a full suspension system to handle varied terrain with ease.

This scooter combines durability with luxury, offering a deluxe high-back seat, full lighting package, and advanced safety features. Whether for daily use or longer recreational outings, the Pathrider 140 XL delivers a smooth, stable ride and the confidence to travel independently.

Key Features

- Heavy-duty 4-wheel scooter for superior outdoor performance
- Full suspension and large pneumatic tyres for a comfortable ride
- Powerful motor with long-range capability
- High-back deluxe seat with headrest and adjustable armrests
- Full LED lighting system, indicators, and rearview mirrors for safety

- Delta tiller with ergonomic controls for ease of use
- Ideal for long-distance travel, hilly terrain, and varied outdoor surfaces
- Supports users seeking a dependable, comfortable, and high-capacity scooter

Specifications

Range Per Charge	Up to 35 km
Maximum Weight Capacity	181 kg
Turning Radius	210 cm right 208 cm left
Maximum Speed	10 km/h
Length & Width	L: 142.8 cm W: 69.2 cm
Seat Type Depth & Width (Standard)	Pillow deluxe top, high back seat W: 49.5 cm (19.5") D: 44.5 cm (17.5")
Seat-To-Deck-Height	47 cm
Seat-To-Ground-Height	66.7 cm – 70.5 cm
Ground Clearance	13.3 cm at midframe
Overall Weight (Not INC Battery)	117.7 kg
Weight Of Heaviest Piece	93 kg (main frame)
Suspension	Front & Rear
3W/4W Option	4W
Front Tyres	13" pneumatic
Rear Tyres	13" pneumatic
Battery Requirements	(2) 70AH
Battery Charger	8A, off-board
Battery Weight	23.45 kg each
Colour Options	Red, Silver, yellow

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au