



PRIDE PATHRIDER 140 XL SCOOTER

\$6,828.59 Original price was: \$6,828.59.<u>\$6,146.00</u>Current price is: \$6,146.00.



SKU: PR140XL-4 RED-H

PRODUCT DESCRIPTION

Purpose & Benefits of Pride Pathrider 140 XL Scooter

The **Pride Pathrider 140 XL Scooter** is a heavy-duty, outdoor-ready mobility scooter built for high performance, comfort, and reliability. Designed for users who require long-distance travel and rugged outdoor capability, the Pathrider 140 XL features a powerful motor, large pneumatic tyres, and a full suspension system to handle varied terrain with ease.

This scooter combines durability with luxury, offering a deluxe high-back seat, full lighting package, and advanced safety features. Whether for daily use or longer recreational outings, the Pathrider 140 XL delivers a smooth, stable ride and the confidence to travel independently.

Key Features

- Heavy-duty 4-wheel scooter for superior outdoor performance
- Full suspension and large pneumatic tyres for a comfortable ride
- Powerful motor with long-range capability
- High-back deluxe seat with headrest and adjustable armrests
- Full LED lighting system, indicators, and rearview mirrors for safety



- Delta tiller with ergonomic controls for ease of use
- Ideal for long-distance travel, hilly terrain, and varied outdoor surfaces
- Supports users seeking a dependable, comfortable, and high-capacity scooter

Specifications

Up to 35 km **Range Per Charge**

Maximum Weight Capacity 181 kg

> **Turning Radius** 210 cm right 208 cm left

Maximum Speed 10 km/h

Length & Width L: 142.8 cm W: 69.2 cm

Pillow deluxe top, high back seat W: 49.5 cm (19.5") **Seat Type Depth & Width (Standard)**

117.7 kg

D: 44.5 cm (17.5")

Seat-To-Deck-Height 47 cm

Seat-To-Ground-Height 66.7 cm - 70.5 cm **Ground Clearance** 13.3 cm at midframe

Overall Weight (Not INC Battery)

Weight Of Heaviest Piece 93 kg (main frame)

> Front & Rear Suspension

3W/4W Option 4W

Front Tyres 13" pneumatic **Rear Tyres** 13" pneumatic

Battery Requirements (2) 70AH

Battery Charger 8A, off-board **Battery Weight** 23.45 kg each

Colour Options Red, Silver, yellow