



### **PRIDE PATHRIDER 140 XL SCOOTER**

**~~\$6,829.00~~ Original price was:**

**~~\$6,829.00.~~ ~~\$6,145.73~~ Current price is: \$6,145.73.**



**SKU: PR140XL-4 RED-H**

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits of Pride Pathrider 140 XL Scooter**

The **Pride Pathrider 140 XL Scooter** is a heavy-duty, outdoor-ready mobility scooter built for high performance, comfort, and reliability. Designed for users who require long-distance travel and rugged outdoor capability, the Pathrider 140 XL features a powerful motor, large pneumatic tyres, and a full suspension system to handle varied terrain with ease.

This scooter combines durability with luxury, offering a deluxe high-back seat, full lighting package, and advanced safety features. Whether for daily use or longer recreational outings, the Pathrider 140 XL delivers a smooth, stable ride and the confidence to travel independently.

### **Key Features**

- Heavy-duty 4-wheel scooter for superior outdoor performance
- Full suspension and large pneumatic tyres for a comfortable ride
- Powerful motor with long-range capability
- High-back deluxe seat with headrest and adjustable armrests
- Full LED lighting system, indicators, and rearview mirrors for safety

- Delta tiller with ergonomic controls for ease of use
- Ideal for long-distance travel, hilly terrain, and varied outdoor surfaces
- Supports users seeking a dependable, comfortable, and high-capacity scooter

## Specifications

<b>Range Per Charge</b>	Up to 35 km
<b>Maximum Weight Capacity</b>	181 kg
<b>Turning Radius</b>	210 cm right 208 cm left
<b>Maximum Speed</b>	10 km/h
<b>Length &amp; Width</b>	L: 142.8 cm W: 69.2 cm
<b>Seat Type Depth &amp; Width (Standard)</b>	Pillow deluxe top, high back seat W: 49.5 cm (19.5") D: 44.5 cm (17.5")
<b>Seat-To-Deck-Height</b>	47 cm
<b>Seat-To-Ground-Height</b>	66.7 cm – 70.5 cm
<b>Ground Clearance</b>	13.3 cm at midframe
<b>Overall Weight (Not INC Battery)</b>	117.7 kg
<b>Weight Of Heaviest Piece</b>	93 kg (main frame)
<b>Suspension</b>	Front & Rear
<b>3W/4W Option</b>	4W
<b>Front Tyres</b>	13" pneumatic
<b>Rear Tyres</b>	13" pneumatic
<b>Battery Requirements</b>	(2) 70AH
<b>Battery Charger</b>	8A, off-board
<b>Battery Weight</b>	23.45 kg each
<b>Colour Options</b>	Red, Silver, yellow