




## PRIDE PATHRIDER 140 XLE MOBILITY SCOOTER

**\$9,585.00**



**SKU:** SCT670504

## PRODUCT DESCRIPTION

### Purpose & Benefits

The Pride Pathrider 140 XL is a robust, full-sized mobility scooter designed for powerful outdoor performance and all-day comfort. Ideal for users who need stability, range, and durability, it features a high-torque motor, full suspension, and large pneumatic tyres to tackle challenging terrain with confidence.

With a sleek modern design and luxury seating, the Pathrider 140 XL combines practicality with comfort. Whether navigating hills, uneven surfaces, or long outdoor journeys, this scooter delivers smooth control and high-end performance—making it an ideal choice for active, independent users.

### Features

- Full suspension system for smooth outdoor riding
- Large pneumatic tyres for enhanced ground clearance
- Powerful drivetrain with 4-pole motor
- Rearview mirrors and full lighting package for safety
- Adjustable high-back seat with headrest and armrests

- Delta tiller with ergonomic wraparound handles
- High weight capacity and extended range battery setup
- Ideal for long-distance outdoor use

## **Specifications**

- Maximum Speed: Up to 10 km/h
- Travel Range: Up to 35 km per charge (depending on conditions)
- Motor: 24-volt DC motor, 4-pole
- Battery: Two 12V, deep cycle (55Ah or 75Ah option)
- Tyres: 33 cm (13") pneumatic (front and rear)
- Ground Clearance: 9.5 cm
- Turning Radius: 170.2 cm
- Maximum User Weight: 181 kg
- Seat Type: Reclining high-back with sliders, headrest, and flip-up armrests
- Suspension: Front and rear independent suspension