



## PROVALE DYSPHAGIA CUP

**\$113.00**



**SKU:** DL2005-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Provale Dysphagia Cup is a specially designed drinking aid for individuals with dysphagia, or difficulty swallowing. This innovative cup is engineered to deliver a measured amount of liquid with each sip, which helps to manage and control the intake for individuals who may otherwise be at risk of choking or aspiration.

### FEATURES

- **Controlled Liquid Delivery:** The cup dispenses a predetermined amount of liquid (often available in different measurements like 5ml or 10ml per sip), which can significantly reduce the risk of aspiration for individuals with dysphagia.
- **Easy-to-Use Design:** The cup is designed with an ergonomic handle and a shape that is easy to hold and control, even for users with limited dexterity or hand strength.
- **Durable and Safe Materials:** Constructed from high-quality, food-grade materials that are both durable and safe for use in healthcare settings.
- **Dishwasher Safe:** Easy to clean and dishwasher safe, making it a practical choice for both healthcare facilities and home use.
- **Transparent or Opaque Options:** Available in either transparent or opaque materials, allowing caregivers and users to monitor liquid intake or provide a more discreet drinking experience, respectively.
- **Wide Base for Stability:** The base of the cup is designed to be wide and stable, reducing the risk of spills.

The Provale Dysphagia Cup is an essential tool for individuals with swallowing difficulties, offering a safe and dignified way to maintain hydration and nutrition. It is highly valued in both clinical settings

and home care, enhancing the quality of life for those living with dysphagia.

## **SPECIFICATIONS**

- Size: 237mL
- Types: 5cc (Blue Top), 10cc (Brown Top)