




PROVALE DYSPHAGIA CUP



SKU: DL2005-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Provale Dysphagia Cup is a specially designed drinking aid for individuals with dysphagia, or difficulty swallowing. This innovative cup is engineered to deliver a measured amount of liquid with each sip, which helps to manage and control the intake for individuals who may otherwise be at risk of choking or aspiration.

FEATURES

- **Controlled Liquid Delivery:** The cup dispenses a predetermined amount of liquid (often available in different measurements like 5ml or 10ml per sip), which can significantly reduce the risk of aspiration for individuals with dysphagia.
- **Easy-to-Use Design:** The cup is designed with an ergonomic handle and a shape that is easy to hold and control, even for users with limited dexterity or hand strength.
- **Durable and Safe Materials:** Constructed from high-quality, food-grade materials that are both durable and safe for use in healthcare settings.
- **Dishwasher Safe:** Easy to clean and dishwasher safe, making it a practical choice for both healthcare facilities and home use.
- **Transparent or Opaque Options:** Available in either transparent or opaque materials, allowing caregivers and users to monitor liquid intake or provide a more discreet drinking experience, respectively.
- **Wide Base for Stability:** The base of the cup is designed to be wide and stable, reducing the risk of spills.

The Provale Dysphagia Cup is an essential tool for individuals with swallowing difficulties, offering a safe and dignified way to maintain hydration and nutrition. It is highly valued in both clinical settings

and home care, enhancing the quality of life for those living with dysphagia.

SPECIFICATIONS

- Size: 237mL
- Types: 5cc (Blue Top), 10cc (Brown Top)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au