



QUEEN COMFORT STANDARD HIGH BACK CHAIR



SKU: AGE600

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Queen Comfort High back supportive chair from Age PrepCare is designed for use in any environment where an adjustable floor to seat height is required: Lounge, dining, bedroom.

The Queen comfort model has the additional support of a high backrest, providing additional back support and comfort.

Available in a variety of colour options in both vinyl and fabric.

FEATURES

- Queen Comfort Wide Healthcare Chair
- High Back
- Rear Castors
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder-coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Depth Adjustment
- 350kg Weight Capacity



• 10-Year Warranty

COLOUR & STYLE

Some colours in the image have been discontinued, select a colour from the available list of options.

DIMENSIONS

Standard

Width	510mm	
Depth	450mm	
Floor to Seat Adjustment	440mm - 525mm	
Back Rest Height	750mm	
Arm Rest Height	200mm	
Weight Capacity	350kg	

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brid	ahton	Showroom	١
------	-------	----------	---

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au

NT



Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au