



## QUEEN COMFORT STANDARD HIGH BACK CHAIR

**\$779.00**



**SKU: AGE600**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Queen Comfort High back supportive chair from Age PrepCare is designed for use in any environment where an adjustable floor to seat height is required: Lounge, dining, bedroom.

The Queen comfort model has the additional support of a high backrest, providing additional back support and comfort.

Available in a variety of colour options in both vinyl and fabric.

### FEATURES

- Queen Comfort Wide Healthcare Chair
- High Back
- Rear Castors
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder-coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Depth Adjustment
- 350kg Weight Capacity

- 10-Year Warranty

## COLOUR & STYLE

Some colours in the image have been discontinued, select a colour from the available list of options.

## DIMENSIONS

### Standard

Width	510mm
Depth	450mm
Floor to Seat Adjustment	440mm - 525mm
Back Rest Height	750mm
Arm Rest Height	200mm
Weight Capacity	350kg