

## QUINTRO HANDYBAR

\$53.00



SKU: 31003



## PRODUCT DESCRIPTION

### QUINTRO HANDYBAR PURPOSE & BENEFITS

The Quintro Handybar Car Aid offers support to help you out of the car safely and securely. The Handy bar inserts into the U-shaped striker plate on your vehicle door frame. The Handy bar 3 in 1 Car Aid uses this strong-point on your vehicle as a simple insertion point that creates an instant, easy to operate "support arm" for exiting the vehicle.

There is a striker on each of the vehicle's door pillars allowing it to be used on both driver and passenger side-doors. The Handy bar provides a leverage point allowing the user to move more easily from a seated position in a vehicle to a standing position with a fraction of the typical difficulty or strain. Located on the Handy bar [handle](#) is a side window breaker.

### FEATURES

- **Promotes Safety & Mobility** - Protects against falls and makes standing or sitting easy.
- **Universal** - Secure fit into passenger or driver side of door latch, no vehicle modification required.
- **No-Slip Ergonomic Grip** - Fits Parallel to car for comfortable grip and optimal leverage.
- **Window Breaker** - Safely and easily breaks glass.
- **Seat Belt Cutter** - Cuts the seat belt in case of an emergency.
- **Weight Capacity:** 136kg
- **Car Types:** Universal fit