



RAIZER II LIFTING CHAIR



SKU: MHLU-105029

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Raizer II Lifting Chair is a mobile lifting chair that helps people who have fallen on the floor to get up again in just a few minutes. Raizer can be operated by a single assistant and requires only limited physical effort beyond the help of a supportive hand.

The Raizer II is easily assembled under a fallen person, and the lifting process is swiftly initiated by a push of a button on the remote control or the control panel on the side of the seat. This means that the assistant does not have to carry out a lift and is exposed to no strain on the back or shoulders as you may experience with other appliances. The helper can with no further assistance handle the Raizer II and help the fallen person.

Our Raizer II Lifting Chair is fitted with safety belts that are to be attached to both backrests to make sure that the fallen person feels secure and comfortable during the lifting process. It is possible to attach additional safety belts on the backrests of Raizer II. Thus, the torso of the fallen person can be secured in several positions for increased comfort.

The Raizer II offers an exclusive and user-friendly design with intuitive sound and light indicators making it easy to quickly assemble the Raizer II at any location. The two backrests are alike, and the four identical legs can be swiftly fitted in any sequence thus making the Raizer II quick and easy to assemble under the fallen person.

FEATURES

- All surfaces of Raizer II are designed to be very easy to clean
- 4 identical legs and 2 identical backrests ensure correct assembly every time

- Intuitive sound and LED indicators for correct assembly
- The seat belt can be removed for washing
- Independently adjustable seat belt placement, on the back rest, ensures improved upper body stability when lifting fallen person
- Purpose-built storage cavity for remote, moulded to the side of seat casing
- Emergency stop button is more intuitive
- Service reminder will indicate when Raizer II needs service
- Battery level indicator
- End stop ensures correct positioning for use and disassembly

OPTIONAL ACCESSORIES

- [Raizer II Headrest](#)

SPECIFICATIONS

Maximum User Weight (SWL):	150kg
Weight	15kg (Seat: 8.5kg / Legs and back supports: 4.5 kg)
Dimensions	68.6cmW x 25.7cmD x 19cmH
Warranty	12 months

VIDEO

<https://www.youtube.com/watch?v=EGmWoly2m3M>

<https://www.youtube.com/watch?v=7Lmr3LPnumE>

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au