



RAIZER II LIFTING CHAIR

~~\$7,299.00~~ Original price was:

~~\$7,299.00.~~\$7,100.00Current price is: \$7,100.00.Sale

Available Online Only!



SKU: MHLU-105029

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Raizer II Lifting Chair is a mobile lifting chair that helps people who have fallen on the floor to get up again in just a few minutes. Raizer can be operated by a single assistant and requires only limited physical effort beyond the help of a supportive hand.

The Raizer II is easily assembled under a fallen person, and the lifting process is swiftly initiated by a push of a button on the remote control or the control panel on the side of the seat. This means that the assistant does not have to carry out a lift and is exposed to no strain on the back or shoulders as you may experience with other appliances. The helper can with no further assistance handle the Raizer II and help the fallen person.

Our Raizer II Lifting Chair is fitted with safety belts that are to be attached to both backrests to make sure that the fallen person feels secure and comfortable during the lifting process. It is possible to attach additional safety belts on the backrests of Raizer II. Thus, the torso of the fallen person can be secured in several positions for increased comfort.

The Raizer II offers an exclusive and user-friendly design with intuitive sound and light indicators making it easy to quickly assemble the Raizer II at any location. The two backrests are alike, and the four identical legs can be swiftly fitted in any sequence thus making the Raizer II quick and easy to assemble under the fallen person.

FEATURES

- All surfaces of Raizer II are designed to be very easy to clean
- 4 identical legs and 2 identical backrests ensure correct assembly every time

- Intuitive sound and LED indicators for correct assembly
- Improved carry bag and trolley
- The seat belt can be removed for washing
- Independently adjustable seat belt placement, on the back rest, ensures improved upper body stability when lifting fallen person
- Purpose-built storage cavity for remote, moulded to the side of seat casing
- Emergency stop button is more intuitive
- Service reminder will indicate when Raizer II needs service
- Battery level indicator
- End stop ensures correct positioning for use and disassembly

OPTIONAL ACCESSORIES

- [Raizer II Headrest](#)

SPECIFICATIONS

Maximum User Weight (SWL):	150kg
Weight	15kg (Seat: 8.5kg / Legs and back supports: 4.5 kg)
Dimensions	68.6cmW x 25.7cmD x 19cmH
Warranty	12 months

VIDEO

<https://www.youtube.com/watch?v=EGmWoly2m3M>

<https://www.youtube.com/watch?v=7Lmr3LPnumE>