



RAIZER M - MANUAL EMERGENCY LIFTING CHAIR

\$3,502.00



SKU: MHLU-107763

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Raizer M is a manual lifting chair that allows a single falls assistant to quickly pick up a fallen person from the floor into a sitting position with a turning handle on the side of the seat.

Ideal for elderly/disabled individuals or simply for those prone to falls at home, the Raizer M is an assistive lifting device designed for floor recovery in domestic and care home environments. The Raizer M Manual Emergency Lifting Chair provides safety and peace of mind for those living in their own home but at risk of falls.

Fall assistance can be safely provided in the comfort of their home without the risk of moving and handling injuries possible from traditional lifting techniques.

FEATURES

- Comfort – The lifting process is very comfortable compared to other solutions on the market.
- Hygiene – Raizer M's materials and surfaces have been designed to be easy to clean and disinfect after use.
- Enable Single-Handed Care – The Raizer M Manual lifting chair can be operated by just one falls assistant, greatly saving on human resources.
- Lifting a person from the floor takes under 5 minutes with the Raizer Manual lift chair.
- Support users up to 150kg

SPECIFICATIONS

Maximum User Weight (SWL):	150kg (23.6 Stone)
Total Lifting Time	3-5 Minutes
Product Weight	12.4 kg

Heaviest Section

7.2 kg

Warranty

2 Years

Dimensions - Standing (L x W x H)

540 x 880 x 1274 mm

VIDEO

[embed]https://www.youtube.com/watch?v=xK7QQJaO_94[/embed]