



RAIZER M - MANUAL EMERGENCY LIFTING CHAIR



SKU: MHLU-107763

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Raizer M is a manual lifting chair that allows a single falls assistant to quickly pick up a fallen person from the floor into a sitting position with a turning handle on the side of the seat.

Ideal for elderly/disabled individuals or simply for those prone to falls at home, the Raizer M is an assistive lifting device designed for floor recovery in domestic and care home environments. The Raizer M Manual Emergency Lifting Chair provides safety and peace of mind for those living in their own home but at risk of falls.

Fall assistance can be safely provided in the comfort of their home without the risk of moving and handling injuries possible from traditional lifting techniques.

FEATURES

- Comfort – The lifting process is very comfortable compared to other solutions on the market.
- Hygiene – Raizer M's materials and surfaces have been designed to be easy to clean and disinfect after use.
- Enable Single-Handed Care – The Raizer M Manual lifting chair can be operated by just one falls assistant, greatly saving on human resources.
- Lifting a person from the floor takes under 5 minutes with the Raizer Manual lift chair.
- Support users up to 150kg

SPECIFICATIONS

Maximum User Weight (SWL):	150kg (23.6 Stone)
Total Lifting Time	3-5 Minutes
Product Weight	12.4 kg

Heaviest Section

7.2 kg

Warranty

12 months

Dimensions - Standing (L x W x H)

540 x 880 x 1274 mm

VIDEO

[embed]https://www.youtube.com/watch?v=xK7QQJaO_94[/embed]

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au