




RAPHAEL X QUATTRO LIFT RECLINE CHAIR



SKU: CHP227640-H

PRODUCT DESCRIPTION

Purpose & Benefits

The Raphael X Quattro lift recline chair is the ultimate comfort and mobility solution, thoughtfully designed to support individuals with limited mobility or postural needs. It empowers users with a smooth sit-to-stand assist and advanced multi-position reclining—ideal for rest, recovery, and day-to-day relaxation.

Powered by the Quattro 4-motor system, users can independently adjust the headrest, backrest, lumbar region, and footrest to achieve precise, ergonomic support. Combined with deep pressure-relieving foam and premium Warwick fabric upholstery, this chair is both a luxurious and functional addition to any living or care environment.

Features

- Advanced Quattro 4-Motor Lift & Recline Mechanism
- Independent control of headrest, backrest, lumbar support, and footrest
- Ergonomic, backlit hand control with memory and home functions
- Deep pressure-relieving foam seat with high-grade cushioning
- Luxurious Warwick Global Suede upholstery

- Quiet and smooth operation for minimal disruption
- Suitable for home and residential aged care use
- Modern aesthetic blends seamlessly with home furnishings
- Back-up battery included for safety during power failure

Specifications

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP227640	Raphael X Quattro Lift Recline Chair - Petite - Golden Grain	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227641	Raphael X Quattro Lift Recline Chair - Petite - Marina Grey	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227642	Raphael X Quattro Lift Recline Chair - Petite - Nordic Stone Vinyl	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227650	Raphael X Quattro Lift Recline Chair - Small - Golden Grain	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm
CHP227651	Raphael X Quattro Lift Recline Chair - Small - Marina Grey	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP227652	Raphael X Quattro Lift Recline Chair - Small - Nordic Stone Vinyl	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm
CHP227660	Raphael X Quattro Lift Recline Chair - Medium - Golden Grain	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227661	Raphael X Quattro Lift Recline Chair - Medium - Marina Grey	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227662	Raphael X Quattro Lift Recline Chair - Medium - Nordic Stone Vinyl	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227670	Raphael X Quattro Lift Recline Chair - Large - Golden Grain	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227671	Raphael X Quattro Lift Recline Chair - Large - Marina Grey	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227672	Raphael X Quattro Lift Recline Chair - Large - Nordic Stone Vinyl	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP227680	Raphael X Quattro Lift Recline Chair - Maxi - Golden Grain	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm
CHP227681	Raphael X Quattro Lift Recline Chair - Maxi - Marina Grey	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm
CHP227682	Raphael X Quattro Lift Recline Chair - Maxi - Nordic Stone Vinyl	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm

*Please note: May require assembly.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au