



## REDGUM 2 IN 1 WALKER/TRANSIT CHAIR



**SKU:** RG4408B-H

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Redgum 2 in 1 Walker/Transit Chair** offers dual functionality as both a seat walker and a transit wheelchair, providing users with flexibility and comfort for indoor and outdoor use. Designed from lightweight aluminium with a durable powder-coated finish, it folds easily for transport and storage, making it ideal for those seeking a versatile and convenient mobility aid.

### Features

- Dual function – converts easily between a seat walker and transit wheelchair
- Durable lightweight aluminium frame with high-quality powder-coated finish
- Side folding design for compact storage and easy transport
- Soft operation hand brakes for smooth control and safety
- Removable and swing-away leg rests for flexibility and comfort
- 115mm wide padded backrest and comfortable 560mm high seat
- Fitted walking stick holder as standard

- 200mm spoked wheels for stability and smooth rolling
- Available in Red, Blue, or White

## Specifications

Specification	Detail
Seat Height	560 mm
Seat Width	450 mm
Seat Depth	310 mm
Overall Width	590 mm
Overall Length	980 mm
Handle Height Range	900 mm - 1000 mm
Unit Net Weight	9.42 kg
Unit Net Weight (No Footrests)	8.3 kg
Wheel Size	200 mm / 8" spoked wheels
Maximum User Capacity	120 kg
Brand	Redgum

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)