



## **REDGUM 2 IN 1 WALKER/TRANSIT CHAIR**

**\$355.00**



**SKU:** RG4408B-H

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The **Redgum 2 in 1 Walker/Transit Chair** offers dual functionality as both a seat walker and a transit wheelchair, providing users with flexibility and comfort for indoor and outdoor use. Designed from lightweight aluminium with a durable powder-coated finish, it folds easily for transport and storage, making it ideal for those seeking a versatile and convenient mobility aid.

### **Features**

- Dual function – converts easily between a seat walker and transit wheelchair
- Durable lightweight aluminium frame with high-quality powder-coated finish
- Side folding design for compact storage and easy transport
- Soft operation hand brakes for smooth control and safety
- Removable and swing-away leg rests for flexibility and comfort
- 115mm wide padded backrest and comfortable 560mm high seat
- Fitted walking stick holder as standard

- 200mm spoked wheels for stability and smooth rolling
- Available in Red, Blue, or White

## Specifications

Specification	Detail
Seat Height	560 mm
Seat Width	450 mm
Seat Depth	310 mm
Overall Width	590 mm
Overall Length	980 mm
Handle Height Range	900 mm - 1000 mm
Unit Net Weight	9.42 kg
Unit Net Weight (No Footrests)	8.3 kg
Wheel Size	200 mm / 8" spoked wheels
Maximum User Capacity	120 kg
Brand	Redgum