



REDGUM 2 IN 1 WALKER/TRANSIT CHAIR

~~\$295.00~~ Original price was: \$295.00. ~~\$265.50~~ Current price is: \$265.50.



SKU: 453146258

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The **REDGUM Dual Purpose Walker / Transit Chair** is ideal for multiple purposes and can be used as either a walker or a transit chair with a small adjustment of the seat rest. The product includes a swing away legrest and walking stick holder.

FEATURES

- Side folding seat walker
- Converts to a transit wheelchair
- Aluminium construction
- Top quality powder coated finish
- Side folding action
- Soft operation hand brake
- Seat height at 56cm from floor
- 11.5cm wide padded backrest
- Removable swing away leg rests fitted
- Walking stock holder fitted as standard
- 20cm spoked wheels

SPECIFICATIONS

- Seat Height - 560mm
- Seat Width - 450mm
- Seat Depth - 310mm
- Unit Overall Width - 590mm

- Unit Overall Length - 980mm
- Handle Height Range - 900mm to 1000mm
- Unit Net Weight - 9.42kg
- Unit Net Weight no footrests - 8.3kg
- Wheel Size - 200mm / 8"
- Maximum User Capacity - 120kg