

## REDGUM ADJUSTABLE BATH BOARD

**\$89.00**



**SKU:** RG9201



### PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

Adjustable showering Bath Board will provide a comfortable position for users who need to be sitting while taking a bath. Adjustable feet -adjust the distance between footrests according to the width of the bathtub, tighten when ready for use. Anti-slip pad on the bottom of the board to improve safety and leaving no damage on the surface of the bathtub.

### BATHBOARD FEATURES

- [Polyurethane](#) self-draining seat
- Removable Backrest that can be fitted on either side
- Rust free construction
- Built-in the soap holder

### SPECIFICATIONS

- Seat width 800mm
- Leg height adjustment range 450mm to 570mm
- Max user capacity 130kg

Bathing can be an uncomfortable and risky activity for those who have mobility issues. Bathers like this Redgum adjustable board help you to feel safe, secure and comfortable while bathing. They are a vital aid for many people with disabilities. We also offer a range of other accessories including bath seats, shower chairs, slipper baths and more!

These amazing Adjustable Bath Boards Provide help to maintain your independence in the comfort of your own home without compromising on quality or safety. With our extensive range of [bathroom](#)

[aids](#) we're sure that we'll have something suitable for all needs - from simple [seat cushions](#) to full-on [wheelchairs](#)!

Buy this Adjustable showering Bath Board online at MediMart using Afterpay or [NDIS method](#) today!