



## REDGUM ADJUSTABLE LEG REST

**\$145.00**



**SKU:** RG63LR-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

This Katie Leg rest is to assist the Disabled or Elderly Finished in Epoxy Coated framework with a Large padded top. The leg rest is height adjustable and can tilt for comfort. Easily stored under the chair when not in use.

### FEATURES

- Can be used as a Leg rest or as a Stool.
- Available in 3 colours.
- Padded platform.
- Height adjustable.
- Angle adjustable.
- Strong lightweight construction.
- Vinyl upholstery.
- Epoxy powder coated metal framework.

### SPECIFICATIONS

- Platform Width: 450mm
- Platform Depth: 380mm
- Platform Height: 325 - 475mm (from floor to top of platform) 25mm adjustment increments
- Product Weight: 5kg
- Maximum User Weight: 120kg (Evenly distributed)