



# **REDGUM ADJUSTABLE LEG REST**

\$145.00 Original price was: \$145.00.<u>\$131.00</u>Current price is: \$131.00.



SKU: RG63LR-H

#### PRODUCT DESCRIPTION

#### **PURPOSE & BENEFITS**

This Katie Leg rest is to assist the Disabled or Elderly Finished in Epoxy Coated framework with a Large padded top. The leg rest is height adjustable and can tilt for comfort. Easily stored under the chair when not in use.

## **FEATURES**

- Can be used as a Leg rest or as a Stool.
- Available in 3 colours.
- Padded platform.
- Height adjustable.
- Angle adjustable.
- Strong lightweight construction.
- Vinyl upholstery.
- Epoxy powder coated metal framework.

## **SPECIFICATIONS**

Platform Width: 450mmPlatform Depth: 380mm

• Platform Height: 325 – 475mm (from floor to top of platform) 25mm adjustment increments

• Product Weight: 5kg

• Maximum User Weight: 120kg (Evenly distributed)