

REDGUM AIR COMFORT CUSHION

~~\$169.00~~ Original price was: ~~\$169.00~~. ~~\$147.00~~ Current price is: **\$147.00**. Sale Available Online Only! While Stocks Last



SKU: RG4040B

PRODUCT DESCRIPTION

Purpose & Benefits

The Redgum Air Comfort Cushion is a pressure-relieving seating solution designed for individuals who require additional comfort and support while seated for extended periods. It uses a combination of foam and an adjustable air cell insert to reduce pressure points and evenly distribute body weight.

This dynamic air and foam hybrid helps improve circulation, reduce the risk of pressure ulcers, and enhance postural support. Ideal for wheelchair users or anyone needing superior sitting comfort, it includes a pump for air level adjustment and a breathable, easy-to-clean cover.

Features

- Hybrid design with foam and air-cell technology
- Adjustable air level for customised support
- Promotes pressure redistribution and improved circulation
- Helps prevent pressure sores and discomfort
- Durable, breathable, and wipe-clean outer cover
- Includes hand pump for easy inflation control

- Lightweight and portable

Specifications

- **Dimensions:** 460mm x 410mm x 100mm
- **Material:** Foam core with internal air cell