



REDGUM ALUMINIUM SOFT BRAKE ADJUSTABLE 8" WALKER - BLUE



SKU: RG4406

PRODUCT DESCRIPTION

Purpose & Benefits

The **Redgum Aluminium Soft Brake Adjustable 8 Inch Walker Blue** is a walking and stability aid designed to support mobility while providing a comfortable seated option when required. It features height adjustable handles and seat, user friendly brakes and park brakes, ergonomic handgrips, a padded seat and backrest, and large 200 mm flat-free tyres for mobility. The lightweight 7 kg design folds for transit or storage and includes a storage bag underneath the seat and a walking stick holder.

Features

- Walking and stability aid.
- Height adjustable seat and handles.
- User friendly brakes and park brakes.
- Ergonomic handgrips.
- Comfortable padded seat.
- Padded backrest.
- Storage bag underneath the seat.

- Walking stick holder.
- 200 mm / 8 inch flat-free tyres.
- 136 kg weight load capacity.
- Lightweight 7 kg design.
- Folds for transit or storage.

Specifications

Overall Dimensions	615 mm (L) x 600 mm (W) x 830 mm (H - backrest)
Overall Folded	370 mm (D) x 600 mm (W) x 870 mm (H)
Seat	405 mm (W) x 290 mm (D) x 500 mm (H)
Wheel Diameter	200 mm / 8 inch
Weight Load Capacity	136 kg
Product Weight	7 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au