




## REDGUM BARIATRIC SEAT WALKER - DARK GREY



**SKU:** RG4413DG

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Redgum Bariatric Seat Walker - Dark Grey** offers strong, stable, and comfortable support for bariatric users. Designed with a wider padded seat, ergonomic hand grips, and large flat-free wheels, it provides a smooth and secure mobility experience both indoors and outdoors.

### Features

- Wider padded seat for maximum comfort
- Strong lightweight construction with 200 kg weight capacity
- Ergonomic hand grips to reduce pressure during use
- Reliable brake and park brake system for safety
- Height adjustable handles and backrest to suit individual users
- Large flat-free wheels for smooth and stable mobility
- Removable storage bag with 5 kg capacity for convenience

- Folds easily for compact storage and transport

## Specifications

- **Overall Dimensions (lowest setting):** 760 mm (W) x 690 mm (L) x 935 mm (H)
- **Seat Dimensions:** 560 mm (W) x 420 mm (L) x 580 mm (H)
- **Folded Dimensions:** 320 mm (W) x 690 mm (L) x 935 mm (H)
- **Bag Size:** 495 mm (W) x 100 mm (L) x 300 mm (H), 5 kg capacity
- **Handle Height Adjustment:** 910 - 1010 mm (5 settings)
- **Wheel Size:** 8" (200 mm) front, 10" (250 mm) rear
- **Product Weight:** 10 kg
- **Maximum User Weight:** 200 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)