



REDGUM BARICARE STEEL QUAD CANE OFFSET HANDLE WITH FOAM GRIP



SKU: RG5611

PRODUCT DESCRIPTION

Purpose & Benefits

The Red Gum BariCare steel quad cane is built for superior support and stability, making it a dependable mobility aid for individuals who need extra assistance while walking. With a heavy-duty steel construction and an offset foam-padded handle, it delivers both strength and comfort in every step.

Its broad quad base enhances balance across a range of surfaces, while the non-slip rubber tips and adjustable height ensure tailored support for different users. Whether for everyday home use or clinical settings, this cane offers reliable weight-bearing capacity and ease of use.

Features

- Heavy-Duty Steel Quad Cane
- Offset Handle with Comfortable Foam Grip
- Wide Quad Base for Increased Stability
- Slip-Resistant Rubber Tips
- Suitable for Home, Hospital, and Aged Care Use

- Australian Designed and Tested

Specifications

- Black Powder-Coated Finish
- Adjustable Height: 740 mm – 965 mm
- 200 mm x 300 mm base
- Safe Working Load: 250 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au