



REDGUM BASIC BARIATRIC COMMODE CHAIR

\$298.00



SKU: B588

PRODUCT DESCRIPTION

Purpose & Benefits

The **Redgum Basic Bariatric Commode Chair** provides a strong, hygienic and practical toileting solution for users who cannot easily access the bathroom. Designed to support up to 280 kg, it combines durability with comfort through its plastic seat, backrest and armrests. The commode is height adjustable and includes a removable bowl, splashguard and backrest, making it ideal for both home and care facility use. Its lightweight, easy-to-clean design ensures safe and convenient daily use for bariatric individuals.

Features

- Allows toileting away from the bathroom for increased convenience and safety
- Strong design with a 280 kg maximum user weight capacity
- Durable plastic seat, backrest and armrests for comfort and hygiene
- Height adjustable through 8 settings (25 mm increments) for user comfort
- Removable bowl with carry handle and lid for easy emptying



- Includes removable splashguard and backrest for versatility and cleaning ease
- Non-slip rubber tips for stability on various surfaces
- Lightweight and easy to move between rooms
- Smooth, easy-to-clean surfaces maintain hygiene standards

Specifications

Bowl

Specification	Details
Model	Redgum Basic Bariatric Commode Chair
Overall Dimensions	685 mm (W) x 620 mm (D) x 700 mm (H) - lowest height setting
Seat Dimensions	590 mm (W) x 470 mm (D) x 410-575 mm (H), adjustable in 25 mm increments (8 settings)
Backrest Pad	400 mm (W) x 150 mm (H); 300 mm from seat to top of backrest
Armrests	50 mm (W) x 310 mm (L) x 210 mm (H) from seat to top; 585 mm between armrests
Product Weight	8.5 kg
Safe Working Load	280 kg
Materials	Plastic seat, backrest and armrests; non-slip rubber tips

Removable with carry handle and lid Additional Features Removable splashguard and backrest; easy-clean design