




REDGUM BASIC BARIATRIC COMMUNE CHAIR



SKU: B588

PRODUCT DESCRIPTION

Purpose & Benefits

The **Redgum Basic Bariatric Commode Chair** provides a strong, hygienic and practical toileting solution for users who cannot easily access the bathroom. Designed to support up to 280 kg, it combines durability with comfort through its plastic seat, backrest and armrests. The commode is height adjustable and includes a removable bowl, splashguard and backrest, making it ideal for both home and care facility use. Its lightweight, easy-to-clean design ensures safe and convenient daily use for bariatric individuals.

Features

- Allows toileting away from the bathroom for increased convenience and safety
- Strong design with a 280 kg maximum user weight capacity
- Durable plastic seat, backrest and armrests for comfort and hygiene
- Height adjustable through 8 settings (25 mm increments) for user comfort
- Removable bowl with carry handle and lid for easy emptying

- Includes removable splashguard and backrest for versatility and cleaning ease
- Non-slip rubber tips for stability on various surfaces
- Lightweight and easy to move between rooms
- Smooth, easy-to-clean surfaces maintain hygiene standards

Specifications

Specification	Details
Model	Redgum Basic Bariatric Commode Chair
Overall Dimensions	685 mm (W) x 620 mm (D) x 700 mm (H) – lowest height setting
Seat Dimensions	590 mm (W) x 470 mm (D) x 410–575 mm (H), adjustable in 25 mm increments (8 settings)
Backrest Pad	400 mm (W) x 150 mm (H); 300 mm from seat to top of backrest
Armrests	50 mm (W) x 310 mm (L) x 210 mm (H) from seat to top; 585 mm between armrests
Product Weight	8.5 kg
Safe Working Load	280 kg
Materials	Plastic seat, backrest and armrests; non-slip rubber tips
Bowl	Removable with carry handle and lid
Additional Features	Removable splashguard and backrest; easy-clean design

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au