




REDGUM BASIC WALKING FRAME



SKU: B4070WS1

PRODUCT DESCRIPTION

Purpose & Benefits of Redgum Basic Walking Frame

The Redgum Basic Walking Frame is designed to provide essential support and stability for individuals who have difficulty walking due to age, injury, or disabilities. This walking frame helps enhance mobility by offering a sturdy structure to lean on, which assists in maintaining balance and reducing the risk of falls. It's particularly useful for indoor use, allowing users to navigate around their homes more safely and independently.

Features

- **Lightweight Design:** Constructed from lightweight aluminium, making it easy to lift and maneuver, even for users with limited strength.
- **Height Adjustable:** Features simple push-button mechanisms that allow the height of the frame to be adjusted to suit the user's comfort and ensure proper posture.
- **Rubber feet:** Equipped with non-slip rubber feet at the base of each leg, providing excellent grip on indoor surfaces and reducing the risk of slipping.
- **Ergonomic Hand Grips:** Includes comfortable hand grips that ensure a secure hold, helping to reduce hand fatigue and improve control while moving.

Specifications

- **Material:** Made primarily from durable, lightweight aluminium for ease of use and longevity.

- **Dimensions:**

- 595mm (Width)
- 530mm (Depth)
- 825mm – 1000mm (Height) 8 settings at 25mm increments.

- **Safe Working Load:** 136kg

- **Weight:** 2.6kg

The Redgum Basic Walking Frame is an excellent mobility aid for those who require a stable support system to aid walking. Its adjustable, lightweight design makes it accessible for a wide range of users, promoting independence and safety in home environments.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au