




REDGUM COMFORTLITE TRANSIT WHEELCHAIR

\$474.00 - \$512.00 Price range: \$474.00 through \$512.00



SKU: RG312-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Discover a new level of comfort and convenience with the Redgum Comfortlite Transit Wheelchair. Specifically crafted for those who require assistance over short to medium distances, this transport chair ensures every journey feels smooth and secure.

FEATURES

- **Compact and Portable:** With its ability to fit snugly in the boot of a car, the Comfortlite is your perfect travel companion.
- **Lightweight Aluminium Frame:** The sturdy aluminium construction allows for a feather-light weight of just 9.5kg, ensuring ease of mobility without compromising on strength.
- **Optimal Support:** Capable of safely accommodating users up to 100kg, you can journey with confidence.
- **Ultimate Comfort:** Luxuriate in the removable padded foam seating and backing, and let the PU-padded fixed armrests support you. The fixed leg rests with flip-up footplates ensure that every journey is as relaxing as it is safe.
- **Foldable Design:** The backrest is not only comfortable but also foldable, accompanied by handy storage pockets.
- **Safety Enhancements:** With double-hand brakes, the wheelchair guarantees additional safety when being pushed by an attendant. Further safety features include a seat belt, removable calf and side straps, and anti-tip wheels at the rear to prevent any untoward incidents.
- **Smooth Rides:** Navigate through any terrain with ease, thanks to the large 12" rear and 6" front flat-free tyres.
- **Aesthetics:** Available in vibrant green or serene blue, ensuring that style accompanies comfort.

SPECIFICATIONS

- 84cm (L) - Footrests up
- 61cm (W) - 27cm (Folded)
- 97.5cm (H) - 73cm (Folded)
- Seat - 48cm (H) x 46cm (W) x 42cm (D)
- Backrest - 45cm (H) x 44.5cm (W)
- Armrest: 20cm (H)
- Wheels:
- 150mm / 6" (Front)
- 300mm / 12" (Rear)
- Product Weight - 9.5kg
- Max. User Weight - 100kg