




REDGUM COMFORTLITE TRANSIT WHEELCHAIR



SKU: RG312-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Discover a new level of comfort and convenience with the Redgum Comfortlite Transit Wheelchair. Specifically crafted for those who require assistance over short to medium distances, this transport chair ensures every journey feels smooth and secure.

FEATURES

- **Compact and Portable:** With its ability to fit snugly in the boot of a car, the Comfortlite is your perfect travel companion.
- **Lightweight Aluminium Frame:** The sturdy aluminium construction allows for a feather-light weight of just 9.5kg, ensuring ease of mobility without compromising on strength.
- **Optimal Support:** Capable of safely accommodating users up to 100kg, you can journey with confidence.
- **Ultimate Comfort:** Luxuriate in the removable padded foam seating and backing, and let the PU-padded fixed armrests support you. The fixed leg rests with flip-up footplates ensure that every journey is as relaxing as it is safe.
- **Foldable Design:** The backrest is not only comfortable but also foldable, accompanied by handy storage pockets.
- **Safety Enhancements:** With double-hand brakes, the wheelchair guarantees additional safety when being pushed by an attendant. Further safety features include a seat belt, removable calf and side straps, and anti-tip wheels at the rear to prevent any untoward incidents.
- **Smooth Rides:** Navigate through any terrain with ease, thanks to the large 12" rear and 6" front flat-free tyres.
- **Aesthetics:** Available in vibrant green or serene blue, ensuring that style accompanies comfort.

SPECIFICATIONS

- 84cm (L) - Footrests up
- 61cm (W) - 27cm (Folded)
- 97.5cm (H) - 73cm (Folded)
- Seat - 48cm (H) x 46cm (W) x 42cm (D)
- Backrest - 45cm (H) x 44.5cm (W)
- Armrest: 20cm (H)
- Wheels:
- 150mm / 6" (Front)
- 300mm / 12" (Rear)
- Product Weight - 9.5kg
- Max. User Weight - 100kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

