




REDGUM ENDURO SUPER LITE CARBON FIBRE SEAT WALKER

\$716.00



SKU: RG4414

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Enduro Seat Walker from Redgum is an innovative mobility aid designed for individuals who require additional support and stability while walking or standing. Featuring a sturdy and durable construction, this seat walker is built to withstand everyday use and is ideal for indoor and outdoor activities.

The Enduro Seat Walker is designed with a lightweight aluminum frame that is both strong and easy to manoeuvre, making it an ideal choice for those who need a mobility aid that is both reliable and functional. The walker features large, solid wheels that provide excellent traction and stability, allowing you to easily navigate over uneven surfaces and rough terrain.

In addition to its sturdy construction, the Enduro Seat Walker is also designed with comfort in mind. The padded seat and backrest provide a comfortable place to rest when you need a break, while the adjustable handles allow you to customise the height of the walker to your specific needs.

FEATURES

- Carbon fibre frame (5.9kg without bag)
- Height adjustable handles
- Ergonomic handgrips
- User friendly loop brake handles
- Brake + park brakes
- Mesh seating
- Padded Backrest
- Large wheels for excellent mobility

- Removable bag with carry handles
- Folds up easily for transit or storage
- Carry handle

SPECIFICATIONS

- Overall: 700mm L x 610mm W x 910mm H (at lowest setting)
- Overall Folded: 700mm L x 250mm W x 770mm H
- Handle Height: 730 – 900mm (12 settings)
- Seat: 440 mm W x 200mm D x 570mm (floor to seat)
- Wheel Diameter: ø200 mm / 8"
- Product Weight: 5.9kg (without bag)
- Max. User Weight: 130kg
- Max. Bag Loading Weight: 5kg