




REDGUM HEIGHT ADJUSTABLE SEAT WALKER

\$245.00 - \$250.00 Price range: \$245.00 through \$250.00



SKU: RG4209-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Redgum Height Adjustable Walker is perfect for indoor and outdoor use for those with limited mobility. It is easy to transport and operate, and has many safety features.

FEATURES

- Aluminium powder coated frame
- Available in three marble finish colours
- Multi height adjustment handles
- Large soft under seat storage bag
- Comfortable folding backrest
- Handbrakes

WHEEL SIZES

6" Walker

- Height adjustable range 470 - 560mm | Handle height adjustment 700 - 900mm.
- Max user weight 120kg | Product Weight 6.5kg

8" Walker

- Height adjustable range 520 - 580mm | Handle height adjustment 750 - 950mm.

- Max user weight 120kg | Product Weight 7.5kg