


REDGUM LIGHTWEIGHT TRI WALKER



SKU: RG4402-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Step into a world of enhanced mobility and daily convenience with the Redgum Lightweight Tri Walker. This modern and sleek design is not only built for stability and comfort but also offers practicality for both daytime errands and night-time strolls.

FEATURES

- **Sturdy Frame:** Crafted with epoxy powder-coated aluminium, ensuring durability and longevity.
- **Contemporary Design:** Modern aesthetics that blend seamlessly with function.
- **Optimal Mobility:** Large wheels coupled with a swivel front wheel make navigating smoother than ever.
- **Storage Solutions:** Comes equipped with a removable zip lock storage bag featuring a carry handle - perfect for everyday essentials.
- **Customisable Comfort:** Enjoy height-adjustable handles and ergonomic rubber hand grips tailored to your needs.
- **Safety First:** With easy-to-engage/release brake levers, locking handbrakes, and reflective sides/front, you're ensured a secure experience.
- **Additional Perks:** Handy holder for your walking stick ensures you have everything you need at arm's length.
- **Hassle-free Assembly:** Set it up in no time, sans tools!
- **Compact Storage:** Designed to fold effortlessly, making storage and transport a breeze.

SPECIFICATIONS

- Overall Dimensions: Length: 580mm, Width: 675mm, Height: 830mm (lowest setting)
- Handle Heights: 830 - 930mm: 5 settings (floor to top of handle)
 - Hole 1: 830mm
 - Hole 2: 855mm
 - Hole 3: 880mm
 - Hole 4: 905mm
 - Hole 5: 930mm
- Handle Width: 480mm (overall)
- Wheels: 190mm Diameter
- Product Weight: 4.6kg (without bag)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au