

## REDGUM OPAL TRANSIT WHEELCHAIR

**\$495.00**



**SKU: RG311**



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Redgum Opal Walker is a transit wheelchair designed to provide stability and support for individuals with mobility challenges. It features a sturdy, lightweight frame and has swing away leg-rests. The walker also includes four large wheels for smooth manoeuvrability, as well as a padded seat and backrest for comfortable seating. With its stylish design and practical features, the Redgum Opal Walker is a reliable and convenient mobility solution for those in need of extra support.

### FEATURES

- Folding frame / fold down back rest
- Padded arm rests
- Removable swing away leg rests
- Adjustable swing up foot plates
- Attendant and occupant hand brakes
- Attendant propelled push handles

### SPECIFICATIONS

- Overall Dimensions: 55 cm (W) × 91 cm (H) × 98.5 cm (D)
- Folded Dimensions: 25 cm (W) × 73.5 cm (H) × 72 cm (D)

- Seat Dimensions: 46 cm (W) × 50 cm (H) × 40 cm (D)
- Arm Rest Height: 24 cm
- Back Rest Height: 41 cm
- Foot Rest Height (from seat): 41–45 cm
- Wheels: 30 cm (back) – 20 cm (front)
- Product Weight: 9.4 kg
- Max. User Weight: 136 kg