

REDGUM OPAL TRANSIT WHEELCHAIR





SKU: RG311

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Redgum Opal Walker is a transit wheelchair designed to provide stability and support for individuals with mobility challenges. It features a sturdy, lightweight frame and has swing away legrests. The walker also includes four large wheels for smooth manoeuvrability, as well as a padded seat and backrest for comfortable seating. With its stylish design and practical features, the Redgum Opal Walker is a reliable and convenient mobility solution for those in need of extra support.

FEATURES

- Folding frame / fold down back rest
- Padded arm rests
- Removable swing away leg rests
- · Adjustable swing up foot plates
- Attendant and occupant hand brakes
- Attendant propelled push handles

SPECIFICATIONS

• Overall Dimensions: 55 cm (W) × 91 cm (H) × 98.5 cm (D)

• Folded Dimensions: 25 cm (W) × 73.5 cm (H) × 72 cm (D)



• Seat Dimensions: 46 cm (W) \times 50 cm (H) \times 40 cm (D)

Arm Rest Height: 24 cmBack Rest Height: 41 cm

Foot Rest Height (from seat): 41-45 cm
Wheels: 30 cm (back) - 20 cm (front)

Product Weight: 9.4 kgMax. User Weight: 136 kg

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au