



REDGUM OVERBED TABLE

\$385.00



SKU: RGOBCT

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The all new Redgum Overbed-Chair table has a new top design that's vinyl wrapped and has a user friendly platform with rounded top edges that are easy to clean. The positive locking adjustable height of the overbed-chair table can be easily adjusted via a large adjustment lever.

Additionally there are two platforms on this product. The main platform is tiltable with 6 angle settings while the fixed side platform features a fixed guard rail.

To ensure that it can fit around chairs the table's base frame is U-shaped while still having a low clearance base to fit under beds.

FEATURES

- User friendly platforms with rounded top edges
- Easy clean sealed platform surface
- Positive locking adjustable height platform
- Large adjustment lever
- Tiltable main platform with retainer - 6 settings
- Fixed side platform with fixed guard rail
- Practical 'U' shaped base frame to fit around chairs
- Low clearance base to fit under beds
- 2x 60mm plastic twin wheel castors
- Beech wood grain platform finish
- Beige epoxy powder coated metal frame

- Chrome mild steel upright