

**REDGUM OVERCHAIR TABLE WITH FLAT TOP****\$285.00****SKU: RG622****PRODUCT DESCRIPTION****PURPOSE & BENEFITS**

The Redgum Overbed-Chair Table is perfect for anyone who wants the convenience of a table that can be easily adjusted to fit their needs.

Whether you're sitting in bed or in a chair, this table provides a stable platform with a variety of adjustment angles to make it perfect for whatever you're using it for.

With its sturdy base and rounded edges, this table is safe and easy to use.

FEATURES

- Plastic platform with retaining lip
- Easy to clean, sealed platform surface
- Locking adjustable height platform
- Large height adjustment lever
- Practical "U" base fits around chairs
- Low clearance base to fit under beds
- Twin wheel locking castors
- Beige epoxy powder coated metal frame
- Chrome mild steel upright

SPECIFICATIONS

- W: 850mm x D: 425mm x H: 800 – 1100mm
- Base Frame:
 - W: 640mm x H: 100mm
- Castors: 50mm – plastic twin wheel locking castors
- Product weight: 11kg

- Max. Load: 12kg – evenly distributed