



REDGUM PRIMO 8" SEAT WALKER



SKU: RG4212BL-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Redgum Primo 8" Seat Walker

The Redgum Primo 8" Seat Walker is a practical and sturdy mobility aid designed to enhance independence and safety for individuals with limited mobility. Featuring a lightweight frame, comfortable seat, and large wheels, this walker is ideal for indoor and outdoor use. Its user-friendly design and adjustable features make it suitable for a wide range of users, providing reliable support for daily activities.

Features

- **Height Adjustable Handles:** Customisable handle heights allow for an ergonomic fit, promoting comfort and proper posture during use.
- **Large 8" Wheels:** Durable wheels provide smooth navigation over various surfaces, ensuring stability and ease of movement indoors and outdoors.
- **Comfortable Padded Seat:** Includes a cushioned seat with a supportive backrest, offering a comfortable resting place when needed.
- **Under-Seat Storage Bag:** A convenient storage bag beneath the seat provides space for personal belongings, making it practical for daily errands.
- **Folding Design:** The walker folds easily for compact storage and transport, making it perfect for travel or space-saving at home.
- **Lightweight and Durable Frame:** The aluminium construction ensures the walker is both lightweight and robust for long-term use.
- **Lockable Hand Brakes:** Features easy-to-use, lockable brakes for enhanced safety and control on slopes or uneven terrain.

Specifications

- **Weight Capacity:** 136 kg
- **Handle Height Range:** Adjustable from 830 mm to 930 mm
- **Seat Height:** 530 mm
- **Walker Width:** 630 mm
- **Unit Weight:** 6.8 kg
- **Wheel Size:** 8 inches
- **Frame Material:** Lightweight aluminium
- **Cleaning:** Wipe clean with a damp cloth

The Redgum Primo 8" Seat Walker is a versatile and reliable mobility solution, offering users the support and convenience they need to maintain independence and enjoy daily activities.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au