



REDGUM QUAD BASE CANE

\$52.00



SKU: RG5620-BLK

PRODUCT DESCRIPTION

Purpose & Benefits of the Redgum Quad Base Cane

The Redgum Quad Base Cane is designed to provide enhanced stability and support for individuals needing assistance with walking. This cane features a four-point base that allows it to stand independently, offering greater balance and security than traditional single-point canes. It is ideal for users with reduced mobility, balance issues, or those recovering from surgery who require extra support to maintain their independence and safety.

Features

- **Quad Base Design:** The four-legged base provides superior stability and allows the cane to stand on its own, reducing the risk of falling and making it easier to retrieve when not in hand.

- **Adjustable Height:** Equipped with an easy-to-use mechanism that allows the cane's height to be adjusted to suit the user's needs, ensuring optimal ergonomics and comfort.
- **Durable Construction:** Made from high-quality materials that ensure durability and can withstand everyday wear and tear.
- **Comfortable Grip:** Features an ergonomically designed handle that fits comfortably in the hand, reducing fatigue and providing a secure grip.
- **Non-slip Rubber Tips:** Each leg is capped with a non-slip rubber tip that enhances grip on various surfaces, providing safe navigation indoors and outdoors.

Specifications

- **Material:** Aluminium, Rubber Tips
- **Safe Working Load:** 100kg

The Redgum Quad Base Cane is an excellent choice for anyone looking for a reliable and stable walking aid. Its quad base design not only offers enhanced stability but also promotes confidence in mobility, making it a favoured option for those with balance concerns or who require additional support while walking.